
Waterpik® Water Flosser Model WP-100

**Hydropropulseur
Waterpik®
Modèle WP-100
Irrigador bucal
Waterpik®
Modelo WP-100**

WELCOME

Traditional Oral Care

Congratulations! By selecting the Waterpik® Water Flosser you are now on your way to better gum health, and a healthier smile!


With a clinically proven combination of water pressure and pulsations, the Waterpik® Water Flosser removes harmful bacteria deep between teeth and below the gumline – where brushing alone won't reach.

We are so convinced that you will be happy with your new Waterpik® Water Flosser, that we back it up with a 14-day guarantee. If you aren't satisfied that your gums are healthier after using the Waterpik® Water Flosser for 14 consecutive days, feel free to return it to us for a full refund.

Thank you for choosing the Waterpik® Water Flosser.


Traditional brushing can't easily reach hidden food debris and harmful bacteria.

Flossing has limited access below the gumline.

Oral Care With The Waterpik® Water Flosser

The Waterpik® Water Flosser removes trapped food debris and harmful bacteria.

Important Safeguards	4
Product Description	6
Getting Started	8
Care of Your Waterpik® Water Flosser	13
Troubleshooting Guide	14
Limited Three-Year Warranty	15

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, the basic safety precautions below should always be followed.

READ ALL INSTRUCTIONS BEFORE USING.

DANGER:

To reduce the risk of electrocution:

- Always unplug product after using.
- Do not handle plug with wet hands.
- Do not immerse in water or other liquid.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub, shower stall or sink.
- Do not reach for a product that has fallen into any liquid. Unplug immediately.
- Check the power supply cord for damage before the first use and during the life of the product.

WARNING:


To reduce the risk of burns, electric shock, fire, serious injury or damage to tissue:



- Do not plug this device into a voltage system that is different from the voltage system specified on the device or charger.
- Do not use this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into any liquid. Contact Water Pik, Inc. or your local distributor which can be found at www.waterpik.com.
- Do not direct water under the tongue, into the ear, nose or other delicate areas. This product is capable of producing pressures that may cause serious damage in these areas. See operating instructions for correct usage.
- Use this product only as indicated in these instructions or as recommended by your dental professional.
- Only use tips and accessories recommended by Water Pik, Inc.
- Fill reservoir with water or other dental professional recommended solution only.
- Do not drop or insert any foreign object into any opening or hose.

- Keep electrical cord away from heated surfaces.
- Do not operate where oxygen or aerosol sprays are being used.
- Do not use mouthwash/rinses that contain iodine, bleach, or tea tree oil.
- Do not use while wearing any oral jewelry. Remove prior to use.
- Do not use if you have an open wound on your tongue or in your mouth.
- If your physician or cardiologist has advised you to receive antibiotic premedication before dental procedures, you should consult your dentist before using this instrument or any other oral hygiene aid.
- For household use only. Do not use outdoors.
- Closely instruct and supervise children and individuals with special needs in the proper use of this product.
- Do not remove the Pik Pocket™ flexible tip from the clear shaft. If it separates for any reason, discard both the tip and shaft and replace with a new Pik Pocket™ Tip. Pik Pocket™ Tip not included in all models.
- Do not use this product for more than 5 minutes in each two-hour period. Longer usage can overheat the motor and cause product failure as well as voiding the warranty.
- Children should be supervised to ensure that they do not play with the appliance.
- This appliance has a polarized plug (one blade is wider than the other). As a safety feature, this plug is intended to fit into a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. Do not attempt to modify the plug in any way.


SAVE THESE INSTRUCTIONS

PRODUCT DESCRIPTION



Plaque
Seeker™
Tip


Classic Jet
Tip


Tongue
Cleaner
Tip


Pik
Pocket™
Tip

Orthodontic
Tip

Toothbrush
Tip

Sonic Tip


NOTE: all tips not included in all models

GETTING STARTED

Plug the cord into an electrical outlet. If the outlet is controlled by a wall switch, be sure it is turned on.


Preparing the Reservoir

Remove reservoir and cover from base. Fill the reservoir with lukewarm water. Replace the reservoir over the base and press down firmly.


Inserting and Removing Tips

Insert tip into the center of the knob at the top of the Water Flosser handle. The colored ring will be flush with the end of the knob if the tip is correctly locked into place. To remove tip from the handle, press the tip eject button and pull tip from handle.

Adjusting the Pressure Setting



Turn the pressure control dial on the base of the unit to the lowest setting (pressure setting No. 1) for first-time use. Gradually increase pressure over time to the setting you prefer – or as instructed by your dental professional.

Placement in Mouth

Lean low over sink and place the tip in mouth. Aim the tip toward teeth.




With your free hand, turn the unit on. Set pressure at desired setting.

Direct the stream by turning the knob at the top of the handle.

Recommended Technique

Direct the stream at a 90-degree angle to your gumline. **Slightly close lips to avoid splashing** but allow water to flow freely from mouth into the sink. For best results, start in the molar area (back teeth) working toward the front teeth. Glide tip along gumline and pause briefly between teeth. Continue until all areas around and between the teeth have been cleaned.


Pause Control

You can temporarily stop the flow any time by pressing the pause button on the handle.

When Finished

Turn the unit off. Empty any liquid left in the reservoir. Unplug the unit.

TIP USAGE

Pik Pocket™ Tip

The Pik Pocket™ Tip is specifically designed to deliver water or anti-bacterial solutions deep into periodontal pockets. It is not intended for general use.

To use the Pik


Pocket™ Tip, lean low over your sink and place the soft tip against a tooth at a 45-degree angle and gently place the tip under the gumline, into the pocket. Turn the unit on to the lowest pressure setting and clean by gently following gumline, inserting tip in between teeth. Continue tracing along the gumline until all areas of the mouth are clean.

NOTE: This tip must be used only with the Pressure Control set to 1. Use on a higher setting will shorten life of tip and power handle

Tongue Cleaner Tip

Place the tongue cleaner on your tongue and start the Water Flosser – use the lowest pressure setting, increase as comfortable.

Gently place the tongue cleaner in the center/midline of your tongue about half way back and with light pressure pull forward.



Orthodontic Tip

The special Orthodontic Tip is uniquely designed to simultaneously brush and rinse hard-to-reach areas around braces and other dental work.

Start on the last molar (back tooth). Aim tip at 90 degree angle to the tooth at the gumline. Gently glide tip along gumline, pausing briefly to lightly brush area between teeth and all around orthodontic bracket, before proceeding to the next tooth.

Continue process until you have cleaned the inside and outside of both the upper and lower teeth.

Toothbrush Tip


Place the toothbrush tip in mouth with the brush head on the tooth at the gum line. The toothbrush tip can be used with or without toothpaste. Turn the Water Flosser on so water flows through the tip. Using a light pressure (bristles should not bend) vibrate the brush back and forth with very short strokes – much like you would with a manual toothbrush. Repeat on adjacent teeth. Water will flow through the brush, to help simultaneously clean and rinse the tooth surface.

Plaque Seeker™ Tip

Direct the jet stream at a 90-degree angle to your gumline. Close lips slightly to avoid splashing but allow water to flow freely from mouth into sink. Start in the molar (back teeth) area working toward the front teeth. Place the tip close to the teeth so that the bristles are gently touching the teeth.

Gently glide tip along gumline, pausing briefly between teeth to gently brush and allow the water to flow between the teeth. Continue until all areas around and between the teeth both inside and outside, top and bottom have been cleaned.

Sonic Tip

Direct the stream at a 90-degree angle to your gumline. **Slightly close lips to avoid splashing** but allow water to flow freely from mouth into the sink. For best results, start in the molar area (back teeth) working toward the front teeth. Glide tip along gumline and pause briefly between teeth for gentle gum massage. Continue until all areas around and between the teeth have been cleaned.

Using Mouthwash and Other Solutions

Your Waterpik® Water Flosser can be used to deliver mouthwash and antibacterial solutions. Ask your dental professional to recommend the appropriate solution for you. Note: use of some antibacterial solutions could shorten the life of your Waterpik® Water Flosser.

After using any special solution in the reservoir, take a moment after every use to rinse the unit. Partially fill the reservoir with warm tap water. Aim the jet tip into your sink. Turn on the unit and let it run until empty. This will help keep the tip from clogging and help your unit to last longer.

Clean the product when necessary by using a soft cloth and mild non-abrasive cleanser.

CARE OF YOUR WATERPIK® WATER FLOSSER

Before cleaning, unplug from the electrical outlet.

Before exposing the Water Flosser to freezing temperatures, remove the reservoir and run the unit until completely empty.

Service Maintenance


Waterpik® Water Flossers have no consumer-serviceable electrical items and do not require routine service maintenance. For all your service needs go to www.waterpik.com. You may also call our toll-free customer service number at 1-800-525-2774 (USA) or 1-888-226-3042 (Canada).

Refer to the serial and model numbers in all correspondence. These numbers are located on the bottom of your unit.

Replacement tips/attachments may be purchased directly from Water Pik, Inc.

Call 1-800-525-2774 (USA) or 1-888-226-3042 (Canada) or visit our web site at www.waterpik.com.

RECOMMENDED TIP REPLACEMENT INTERVALS

TROUBLESHOOTING GUIDE

PROBLEM	CAUSE	CORRECTIVE ACTION
Tips and Tongue Cleaner		
Leakage between the tip and handle knob.	1) Tip not fully engaged.	1) Remove tip and reinstall.
Pik Pocket™ Tip		
Tear in soft tip.	1) Tip is worn.	1) Replace with new tip.
Water Flosser Base		
Hose to the handle leaks.	1) Damage to the hose.	1) Visit www.waterpik.com or contact Water Pik, Inc. customer service for a self service repair kit.
Reservoir leaks.	1) Missing reservoir valve.	1) Replace reservoir valve Contact Water Pik, Inc. customer service for a replacement valve.
	2) Reservoir valve is upside down.	2) Turn valve right-side up.
Inadequate pressure.	1) Reservoir not positioned properly.	1) Re-position reservoir. Press down firmly on the reservoir.
	2) Tip is clogged with debris.	2) Flush tip under faucet or replace tip.
Unit doesn't start.	1) Motor not running.	1) Be sure outlet is functioning; test with another appliance (check electrical reset button and/or wall switch).

Still have questions? Go to www.waterpik.com or call Water Pik, Inc. Customer Service. Toll -Free Number: 1-800-525-2774 (U.S.A.) 1-888-226-3042 (Canada).

LIMITED THREE-YEAR WARRANTY

Water Pik, Inc. warrants to the original purchaser/owner of this new product that it is free from defects in materials and workmanship for three years from date of purchase. Save your receipt as your proof of purchase date. We will replace any part of the product, which in our opinion is defective, provided the product has not been abused, misused, altered or damaged after purchase and was used according to these instructions. This limited warranty excludes accessories or consumable parts such as tips, etc.

This warranty gives you specific legal rights. You may also have other rights that vary depending on the law in your place or residence.

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Water Pik, Inc.

1730 East Prospect Road
Fort Collins, CO 80553-0001 USA

www.waterpik.com

Waterpik® is a trademark of Water Pik, Inc. registered in Argentina, Australia, Canada, Chile, China, Czech Republic, Hungary, India, Mexico, New Zealand, Pakistan, Poland, Russian Federation, Singapore, South Africa, Switzerland, Taiwan, Ukraine, United Kingdom, and the United States.

Waterpik (stylized) is a trademark of Water Pik, Inc. registered in Canada, EU, Japan, Mexico, and the United States.

treat yourself better™, Pik Pocket™ and Plaque Seeker™ are trademarks of Water Pik, Inc.

Form No. 20013283-F AA

©2009 Water Pik, Inc.

Designed in USA.

Conçu aux É.-U.

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Made in CHINA.

Fabriqu   en CHINE.

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JOB INFORMATION

Product Description: WaterPik

Instruction Manual Model WP-100

Client: WaterPik

Manufacturer:

Blueprint Ref: 1007 file, 5x7 folded

UPC No:

Last Revision Date: 1/05/10

STUDIO INFORMATION

Job No: 2272.0007.005

File No: 2272.0007.005 WP 100 IG.indd

Studio Contact: Client Service: Mark Morse

Creative Lead: Scott Smith

Production Mgr: Doug Bennett

CLIENTS PACKAGING APPROVALS

INITIALS	DATE
BRAND	
LEGAL	

IMPORTANT NOTICE

brandimage has checked this artwork for accuracy. Final approval is the client's responsibility. Please double check for errors before production. We recommend that color proofs be submitted to us for approval before a production.

TECHNICAL INFORMATION

Application/Version: Adobe Indesign CS3

Layers/Pages: 5/1

Pages/Images:

Fonts Used: Refer to file

Signatures:

Name

Release Date: / /

Approved For Release

COMPUTER EQUIVALENTS

TINTS / BUILT COLORS

brandimage APPROVALS

gmb	/ /	/ /
Art/Operator	Date	Production Director
Account Executive	Date	Proofreader
Design Director	Date	Proofreader